



Product Name	Bottle-Top Dispenser	New Bottle-Top Dispenser	Electronic Bottle-Top Dispenser	Electronic Titrator	Vacuum Aspiration System	Economical Vacuum Aspirator
Model	DispensMate	DispensMate-Pro/S	dFlow	dTrite	SAFEVAC	EcoVAC
Application	Reagent pipetting and liquid dispensing			Used in titration operations in chemical analysis, food industry, environmental analysis	Liquid collection and storage treatment, including cell culture medium, suspension, supernatant.	
Operation	Manual	Manual	electric	electric	electric	electric
Volume Range	0.5-5ml 1.0-10ml 2.5-25ml 5.0-50ml	0.5-5ml 1.0-10ml 2.5-25ml 5.0-50ml 10-100ml	0.1ml-99.9ml	0.01ml-99.99ml	Vacuum range: 0-600mbar Bottle volume:4L	Vacuum range: 0-500mbar Bottle volume:2L
Reagent Recovery Function	No	Yes	Yes	Yes	No	No
Features		-S Enhanced seal -Pro Glass piston	-	-	Sensitive level sensor,Self-closing connector	1L bottle optional, Switch manual and continuous aspiration
Fully autoclavable	The part where the liquid flows	Whole	-	-	The part where the liquid flows	The part where the liquid flows
Accessories	Screw adapter	Screw adapter	Spiral dispensing pipe	Spiral dispensing pipe	Handle adapter	Handle adapter