

Series	Image	Model	Display	Max. speed(RCF)	Rotor Capacity
palm micro centrifuges		D1008	-	7000rpm (2680×g)	8×2ml 32×0.2ml
High speed Mini Centrifuge		D2012Plus	LCD	15000rpm 15100×g)	12×0.2/0.5 /1.5/2ml
High speed Micro Centrifuge		D3024	LCD	15000rpm (21380×g)	24×1.5/2ml, 36×0.5ml, 4×PCR eight strips, 12×5ml(cone), 18×5ml(round),
High-speed Refrigerated Micro Centrifuge		D1524R	LCD	15000rpm (21380×g)	
hematocrit centrifuge		DM1424	LCD	14000 rpm (18620×g)	24×1.5/2ml(14000rpm), 24×capillary tubes (12000 rpm)
		DM1224	LCD	12000 rpm (13680×g)	24×capillary tubes
DM0506 Series		DM0506	LCD	5000rpm (2600×g)	6×15ml/10ml /7ml/1.5-5ml
		DM0412	LCD	4500rpm (2490×g)	8×15ml 12×10ml/7ml 1.5-5ml
		DM0412P (PRP/PRF/CGF)	LCD	4500rpm (2490×g)	8×15ml 12×10ml/7ml 1.5-5ml
DM0424 & DM0408 Series		DM0408	LED	4000rpm (2500×g)	8×15ml 12×10ml/7ml 1.5-5ml
		DM0424	LCD	4000rpm (2500×g)	4×50ml 24×15ml
DM0636 Series		DM0636	LCD	6000rpm (4020×g)	4×100ml 36×3-7ml