











Series	Image	Model	Application	Operation Method	Volume Range	Channel Option	Autoclavable
Electronic Pipette		dPette+	Aqueous solution	Electronic	0.5-1000ul	Single Channel Adj.	YES (lower part)
		dPette	Aqueous solution	Electronic	0.5-1000ul	Single Channel Adj.	YES (lower part)
		MicroPette	Aqueous solution	Manual	0.1ul-10ml	Single Channel Adj. & Fixed; 8 Channel Adj.; 12 Channel Adj.;	YES (lower part)
		MicroPette Plus	Aqueous solution	Manual	0.1ul-5000ul	Single Channel Adj. & Fixed; 8 Channel Adj.; 12 Channel Adj.;	YES
		MiniPette	Aqueous solution	Manual	25-200ul	Single Channel Fixed	NO
Pipette Filler		Levo Plus	Aqueous solution	Electronic	0.1-100ml	Single channel	YES (silicon adapter&nozzle)
		Levo Me	Aqueous solution	Electronic	0.1-100ml	Single channel	YES (silicon adapter&nozzle)
Pipette Controller		Levo	Aqueous solution	Manual	0.1-100ml	Single channel	YES (silicon adapter&nozzle)
Pipette Pump		Levo E	Aqueous solution	Manual	2-25ml	Single channel	NO
Stepper		StepMate	Aqueous solution, high viscosity or volatile solution	Manual	0.5-50ml	Single channel	NO