

Image					
Models	OS20-Pro	OS20-S	OS40-Pro	OS40-S	OS70-Pro
Max. Stirring Quantity[H ₂ O]	20L	20L	40L	40L	70L
Constant Speed	YES	NO	YES	NO	YES
Display	LCD	LED	LCD	LED	LCD
Speed Range	50-2200rpm	50-2200rpm	50-2200rpm	50-2200rpm	50-1100rpm
Max. Torque	40N•cm	40N•cm	60N•cm	60N•cm	300N•cm
Viscosity Max.	10000mPa•s	10000mPa•s	50000mPa•s	50000mPa•s	100000mPa•s