

Model	Image	Function	Temp. Setting Range	Mixing speed	Blocks
HB60-S		Heating	25°C,30°C, 37°C,45°C,60°C,	-	5ml×30+15ml×8
HB120-S		Heating	25°C-120°C	-	0.2ml×54 0.5ml×40 1.5/2ml×40 5/15ml×28 5ml blood collection tube×45 50ml×8 (0.2ml+0.5ml+1.5/2ml)×18 0.2ml×96 Microplate Bead/water bath
HB105-S1/ S2		Heating	25°C-105°C	-	0.2ml×30 0.5ml×20 1.5ml×20
HB150-S1/ S2		Heating	25°C-150°C	-	2ml×20 15ml×12 50ml×4 Bead/water bath
H100-Pro		Heating	25°C-100°C	-	
HC110-Pro		Heating& Cooling	-5°C-110°C	-	0.5ml×24 1.5ml×24 2ml×24 5ml×5 15ml×8
HCM100-Pro		Heating& Cooling& Mixing	0°C-100°C	200- 1500rpm	50ml×4 0.2ml×96 Microplate
HM100-Pro		Heating& Cooling	15°C-100°C	200- 1500rpm	
Mini H100		Heating	25°C-100°C	-	0.2ml×40 0.5ml×24
Mini HC100		Heating& Cooling	0°C-100°C	-	1.5ml×15 2ml×15 5ml×6 15ml×4 50ml×2
Mini HCL100		Heating& Cooling	0°C-100°C	-	