

hemochroma PLUS

The simplest way of testing hemoglobin



hemochroma^{PLUS}

- Results in 3 seconds
- Cyanide-free cuvette
- Capillary or tube blood
- Battery or power from USB
- Test instructions in color LCD
- Long open-stability

Boditech Med Inc.

43, Geodudanji 1-gil, Dongnae-myeon,
Chuncheon-si, Gang-won-do 24398, Korea
Tel: +82-33-243-1400, Fax: +82-33-243-9373
sales@boditech.co.kr / www.boditech.co.kr



Find out more on www.boditech.co.kr



hemochroma^α PLUS


Simple and easy way testing hemoglobin



> Specifications

Specification	Value
Principle	Spectrophotometry
Sample	Whole blood
Sample volume	12 – 15 uL
Detection range	0.5 – 27 g/dL
Reading time	Approx. 3 sec.
Data storage	1,000 results
Power	4 x AA batteries or power from USB
Dimensions	106 x 151 x 39 mm
Weight	215 g
Calibration	Calibration information in ID chip

> Performance



> Applications

- Screening anemia
- Blood banks
- Transfusion centers
- Mother and child health centers
- Dialysis centers

> Simple procedure


Sample collection



Measurement



or



> Ordering information

Product	Cat#	Unit
hemochroma PLUS	FPRR005	1 set
Microcuvette	CFPO-27	50 ea/bottle, 4 bottles/box