

Series	Model	Image	Stirring Position	Display	Work plate Size	Max. Heating Temp.	Max. stirring quantity[H <sub>2</sub> O]
550 series	MS7-H550-Pro		1	LCD	184x184mm	550°C	20L
	MS7-H550-S		1	LED	184x184mm	550°C	10L
380 series	MS-H380-Pro		1	LCD	140x140mm	380°C	20L
340 series	MS-H-Pro <sup>T</sup> (with timer)		1	LCD	φ135mm	340°C	20L
	MS-H-Pro <sup>A</sup>		1	LCD	φ135mm	340°C	20L
	MS-H-Pro <sup>+</sup>		1	LCD	φ135mm	340°C	20L
	MS-H-S		1	Scale	φ135mm	340°C	20L
	MS-H340-S		4	LCD	φ135mm	340°C	10L for single position
280 series	MS-H280-Pro		1	LED	φ135mm	280°C	3L
120 series	MS-H-S10		10	Scale	180x450mm	120°C	0.4Lx10